

DíamoDent[®]

RESTORATIVE CATALOG AND PRICE LIST



MARCH 2013
NINTH EDITION

Table of Contents

Prosthetic products compatible with:

Nobel Biocare®

Internal Connection - NobelReplace® Tapered & Straight Groovy®, NobelSpeedy®

Replace, Replace® Select Tapered & Straight Systems

3.5mm Narrow Platform (NP) implants®	5-6
4.3mm Regular Platform (RP) implants®	7-8
5.0mm Wide Platform (WP) implants®	9-10
6.0mm implants®	11-12

External Connection – Branemark® System

3.75 & 4mm Regular Platform (RP) implants®	43-44
5mm Wide Platform (WP) implant®	48

Zimmer®

3.3mm & 3.7mm (3.5mm Platform) Screw-Vent Zimmer implants®	13-14
4.7mm (4.5mm Platform) Screw-Vent Zimmer implant®.....	15-16
6mm (5.7mm Platform) Screw-Vent Zimmer implant®	17-18

Table of Contents

Implant Direct®

RePlant® Implant System

3.5mm Diameter (3.5mm Platform) implant	5-6
4.3mm Diameter (4.3mm Platform) implant	7-8
5.0mm Diameter (5.0mm Platform) implant	9-10
6.0mm Diameter (6.0mm Platform) implant	11-12

RePlus® and ReActive® Implant Systems

3.7mm Diameter (3.5mm Platform) implants	5-6
4.7mm Diameter (4.3mm Platform) implants	7-8
5.7mm Diameter (5.0mm Platform) implants	9-10

Legacy® Implant System

3.7mm Diameter (3.5mm Platform) implant	13-14
4.7mm Diameter (4.5mm Platform) implant	15-16
5.7mm Diameter (5.7mm Platform) implant	17-18

BioHorizons®

Internal Hex Connection

3.5mm (3.5mm Platform) implants®	19-20
4.0mm (4.5mm Platform) implant®.....	21-22
5.0mm & 6.0mm (5.7mm Platform) implant®	23-24

Table of Contents

BIOMET® 3i®

Internal Hex Connection – 3i Certain Implants

3.4mm Platform Certain implants®	25-26
4.1mm Platform Certain implants®	27-28
5.0mm Platform Certain implants®	29-30
6.0mm Platform Certain implants®	31-32

External Hex Connection

3.25 & 3.3mm MicroMiniplant implants®	42
Osseotite XP 3/4, 3.25 & 3.3mm Miniplants, 3.75 & 4mm Standard implants®	43-44
5mm implants	45-46
6mm implants	47

Lifecore® / Keystone Dental®

External Hex Connection – RESTORE® Implant System

3.3 & 3.4mm Small Diameter (SD) implants®	49
3.75, 4 & 4.2mm Regular Diameter (RD) implants®	43-44
5 & 6mm Wide Diameter (WD) implants®	50

Straumann®

4.1mm/4.8mm Regular Neck (RN) implants®	33-34
4.8mm Wide Neck (WN) implants®	35-36
3.3mm Narrow Neck (NN) implants®	37

Table of Contents

Friadent®

3.4mm Frialit-2 & XiVE implants®	38
3.8 Frialit-2 & XiVE implants®	39
4.5mm Frialit-2 & XiVE implants®	40
5.5mm Frialit-2 & XiVE implants®	41

Steri-Oss®/ Nobel Biocare®

3.8 & 4.5mm HL implants®	43-44
4.3mm Replace implant®	51
5mm Replace and HL Threaded implants®	52

Tools and Accessories

Screwdrivers	53
Hader Bars	53
Driver Tips	54
Lab Screwdrivers	55

Price List	56-60
-------------------------	-------

Patents

Some of the products in this catalog are covered by the following patents:

US 7,104,797 B2
US 6,923,648 B1
US D562,453 S

DiamoDent®

Compatible with: -3.5mm Narrow Platform (NP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® SelectTapered & Straight implants -3.5mm Diameter (3.5mm Platform) RePlant® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Purple

AA100 Engaging, 0.5mm Collar, 3.7mm Profile
AA101 Non-Engaging, 0.5mm Collar, 3.7mm Profile

Requires ComScrew catalog number AA660



AA100



AA101

UCLA ABUTMENT, Clear Plastic, Purple

AA200 Engaging, 1.0mm Collar, 4.3mm Profile
AA201 Non-Engaging, 1.0mm Collar, 4.3mm Profile

Requires ComScrew catalog number AA660



AA200



AA201

OVERDENTURE ABUTMENT, Clear Plastic, Purple

AA401 Non-Engaging

Requires ComScrew catalog number AA660

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

AA505 1mm Collar, 4.5mm Profile
AA510 2mm Collar, 4.5mm Profile
AA515 4mm Collar, 4.5mm Profile
AA520 2mm Collar, 5.3mm Profile
AA525 4mm Collar, 5.3mm Profile
AA535 3mm Collar, 5.3mm Profile

Requires ComScrew catalog number AA660



AA505



AA510



AA515



AA520



AA525



AA535

ComScrew®, Titanium

AA660 UniGrip®
35 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: -3.5mm Narrow Platform (NP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® SelectTapered & Straight implants
-3.5mm Diameter (3.5mm Platform) RePlant® Implant Direct® implants

GUIDE PIN/WAXING SCREW, Titanium, Purple

AA680

UniGrip®



THREADED IMPRESSION COPING, Titanium, Purple

AA700

Closed Tray, 4.1mm Profile



THREADED IMPRESSION COPING, Titanium, Purple

AA750

Open Tray, 4.2mm Profile



TEMPORARY ABUTMENT, Titanium

AA800

Engaging, 1mm Collar, 4.5mm Profile

AA801

Non-Engaging, 1mm Collar, 4.5mm Profile

Requires ComScrew catalog number AA660



AA800



AA801

IMPLANT REPLICA, Aluminum, Purple

AA910



ZIRCONIA ABUTMENT BASE, Titanium, Purple

AA1000

Engaging, 0.35 Collar, 4.0mm Profile

Requires ComScrew catalog number AA660



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 4.3mm Regular Platform (RP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 4.3mm Diameter (4.3mm Platform) RePlant® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Yellow

AB100 Engaging, 0.5mm Collar, 4.5mm Profile
AB101 Non-Engaging, 0.5mm Collar, 4.5mm Profile

Requires ComScrew catalog number AB660



AB100



AB101

UCLA ABUTMENT, Clear Plastic, Yellow

AB200 Engaging, 1.0mm Collar, 5.0mm Profile
AB201 Non-Engaging, 1.0mm Collar, 5.0mm Profile

Requires ComScrew catalog number AB660



AB200



AB201

OVERDENTURE ABUTMENT, Clear Plastic, Yellow

AB401 Non-Engaging

Requires ComScrew catalog number AB660

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

AB505 1mm Collar, 5.3mm Profile
AB510 2mm Collar, 5.3mm Profile
AB515 4mm Collar, 5.3mm Profile
AB520 2mm Collar, 6.1mm Profile
AB525 4mm Collar, 6.1mm Profile
AB535 3mm Collar, 6.1mm Profile

Requires ComScrew catalog number AB660



AB505



AB510



AB515



AB520



AB525



AB535

ComScrew® Titanium

AB660 UniGrip®

35 Ncm



DíamoDent®

Compatible with: - 4.3mm Regular Platform (RP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 4.3mm Diameter (4.3mm Platform) RePlant® Implant Direct® implants

GUIDE PIN/WAXING SCREW, Titanium

AB680

UniGrip®



THREADED IMPRESSION COPING, Titanium, Yellow

AB700

Closed Tray, 4.5mm Profile



THREADED IMPRESSION COPING, Titanium, Yellow

AB750

Open Tray, 4.7mm Profile



TEMPORARY ABUTMENT, Titanium

AB800

Engaging, 1mm Collar, 5.3mm Profile

AB801

Non-Engaging, 1mm Collar, 5.3mm Profile



AB800



AB801

Requires ComScrew catalog number AB660

IMPLANT REPLICA, Aluminum, Yellow

AB910



ZIRCONIA ABUTMENT BASE, Titanium, Yellow

AB1000

Engaging, 0.35 Collar, 4.8mm Profile

Requires ComScrew catalog number AB660



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with:

- 5.0mm Wide Platform (WP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 5.0mm Diameter (5.0mm Platform) RePlant® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Blue

AC100	Engaging, 0.5mm Collar, 5.2mm Profile
AC101	Non-Engaging, 0.5mm Collar, 5.2mm Profile

Requires ComScrew catalog number AB660



AC100



AC101

UCLA ABUTMENT, Clear Plastic, Blue

AC200	Engaging, 1.0mm Collar, 5.5mm Profile
AC201	Non-Engaging, 1.0mm Collar, 5.5mm Profile

Requires ComScrew catalog number AB660



AC200



AC201

OVERDENTURE ABUTMENT, Clear Plastic, Blue

AC401	Non-Engaging
-------	--------------

Requires ComScrew catalog number AB660

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

AC505	1mm Collar, 6.0mm Profile
AC510	2mm Collar, 6.0mm Profile
AC515	4mm Collar, 6.0mm Profile
AC520	2mm Collar, 6.8mm Profile
AC525	4mm Collar, 6.8mm Profile

Requires ComScrew catalog number AB660



AC505



AC510



AC515



AC520



AC525

ComScrew® Titanium

AB660	UniGrip®
-------	----------

35 Ncm

9

CUSTOMER SERVICE
888-281-8850



DíamoDent®

Compatible with:

- 5.0mm Wide Platform (WP) NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 5.0mm Diameter (5.0mm Platform) RePlant® Implant Direct® implants

GUIDE PIN/WAXING SCREW, Titanium

AB680

UniGrip®



THREADED IMPRESSION COPING, Titanium, Blue

AC700

Closed Tray, 5.5mm Profile



THREADED IMPRESSION COPING, Titanium, Blue

AC750

Open Tray, 5.5mm Profile



TEMPORARY ABUTMENT, Titanium

AC800

Engaging, 1mm Collar, 5.4mm Profile

AC801

Non-Engaging, 1mm Collar, 5.4mm Profile



AC800



AC801

Requires ComScrew catalog number AB660

IMPLANT REPLICA, Aluminum, Blue

AC910



ZIRCONIA ABUTMENT BASE, Titanium, Blue

AC1000

Engaging, 0.35 Collar, 5.5mm Profile

Requires ComScrew catalog number AB660



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 6.0mm NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 6.0mm Diameter (6.0mm Platform) RePlant® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Green

AD100	Engaging, 0.5mm Collar, 6.2mm Profile
AD101	Non-Engaging, 0.5mm Collar, 6.2mm Profile

Requires ComScrew catalog number AB660



AD100



AD101

UCLA ABUTMENT, Clear Plastic, Green

AD200	Engaging, 1.0mm Collar, 6.5mm Profile
AD201	Non-Engaging, 1.0mm Collar, 6.5mm Profile

Requires ComScrew catalog number AB660



AD200



AD201

CEMENT-ON POST, Titanium

AD510	2mm Collar, 7.0mm Profile
-------	---------------------------

Requires ComScrew catalog number AB660



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 6.0mm NobelReplace® Tapered & Straight Groovy®, NobelSpeedy® Replace, Replace® Select Tapered & Straight implants
- 6.0mm Diameter (6.0mm Platform) RePlant® Implant Direct® implants

ComScrew® Titanium

AB660

UniGrip®

35 Ncm



GUIDE PIN/WAXING SCREW, Titanium

AB680

UniGrip®



THREADED IMPRESSION COPING, Titanium, Green

AD700

Closed Tray, 6.5mm Profile



IMPLANT REPLICA, Aluminum, Green

AD910




ZIRCONIA ABUTMENT BASE, Titanium, Green

AD1000

Engaging, 0.35 Collar, 6.5mm Profile

Requires ComScrew catalog number AB660



DíamoDent®

Compatible with: - 3.3mm & 3.7mm (3.5mm Platform) Screw-Vent Zimmer® implants
- 3.7mm Diameter (3.5mm Platform) Legacy® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Green

P110 Hexed, 0.8mm Collar, 4.5mm, Profile, TiteFit™
P101 Non-Hexed, 0.8mm Collar, 4.5mm Profile

Requires ComScrew catalog number P655



P110



P101

UCLA ABUTMENT, Clear Plastic Green

P200 Hexed, 0.9mm Collar, 4.3mm Profile
P201 Non-Hexed, 0.9mm Collar, 4.3mm Profile

Requires ComScrew catalog number P655



P200



P201

OVERDENTURE ABUTMENT, Clear Plastic, Green

P401 Non-Hexed

Requires ComScrew catalog number P655

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

P505 1mm Collar, 4.5mm Profile, TiteFit™
P510 2mm Collar, 4.5mm Profile, TiteFit™
P515 4mm Collar, 4.5mm Profile, TiteFit™
P520 2mm Collar, 5.5mm Profile, TiteFit™
P525 4mm Collar, 5.5mm Profile, TiteFit™
P535 3mm Collar, 5.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



P505



P510



P515



P520



P525



P535

ComScrew®, Titanium

P655 .050" Hex

30 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 3.3mm & 3.7mm (3.5mm Platform) Screw-Vent Zimmer® implants
- 3.7mm Diameter (3.5mm Platform) Legacy® Implant Direct® implants

GUIDE PIN/WAXING SCREW, Titanium

P680 .050" Hex



THREADED IMPRESSION COPING, Titanium, Green

P700 Closed Tray, 4.6mm Profile



THREADED IMPRESSION COPING, Titanium, Green

P750 Open Tray, 4.7mm Profile



TEMPORARY ABUTMENT, Titanium

P800 Hexed, 1mm Collar, 4.5mm Profile
P801 Non-Hexed, 1mm Collar, 4.5mm Profile



P800

P801

Requires ComScrew catalog number P655

IMPLANT REPLICA, Stainless Steel

P900

Suitable for soldering



IMPLANT REPLICA, Aluminum, Green

P910




ZIRCONIA ABUTMENT BASE, Titanium, Green

P1000 Hexed, 0.35 Collar, 4.0mm Profile
Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 4.7mm (4.5mm Platform) Screw-Vent Zimmer® implants
- 4.7mm Diameter (4.5mm Platform) Legacy® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Purple

R110 Hexed, 0.8mm Collar, 4.7mm Profile, TiteFit™
R101 Non-Hexed, 0.8mm Collar, 4.7mm Profile

Requires ComScrew catalog number P655



R110



R101

UCLA ABUTMENT, Clear Plastic, Purple

R200 Hexed, 0.7mm Collar, 4.6mm Profile
R201 Non-Hexed, 0.7mm Collar, 4.6mm Profile

Requires ComScrew catalog number P655



R200



R201

OVERDENTURE ABUTMENT, Clear Plastic, Purple

R401 Non-Hexed

Requires ComScrew catalog number P655

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

R505 1mm Collar, 5.5mm Profile, TiteFit™
R510 2mm Collar, 5.5mm Profile, TiteFit™
R515 4mm Collar, 5.5mm Profile, TiteFit™
R520 2mm Collar, 6.5mm Profile, TiteFit™
R525 4mm Collar, 6.5mm Profile, TiteFit™
R535 3mm Collar, 6.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



R505



R510



R515



R520



R525



R535

ComScrew™, Titanium

P655 .050" Hex
30 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: - 4.7mm (4.5mm Platform) Screw-Vent Zimmer® implants
- 4.7mm Diameter (4.5mm Platform) Legacy® Implant Direct® implants

GUIDE PIN/WAXING SCREW, Titanium

P680 .050" Hex



THREADED IMPRESSION COPING, Titanium, Purple

R700 Closed Tray, 4.9mm Profile



THREADED IMPRESSION COPING, Titanium, Purple

R750 Open Tray, 5.2mm Profile



TEMPORARY ABUTMENT, Titanium

R800 Hexed, 1mm Collar, 5.4mm Profile
R801 Non-Hexed, 1mm Collar, 5.4mm Profile



R800

R801

IMPLANT REPLICA, Stainless Steel

R900
Suitable for soldering



IMPLANT REPLICA, Aluminum, Purple

R910




ZIRCONIA ABUTMENT BASE, Titanium, Purple

R1000 Hexed, 0.35 Collar, 5.0mm Profile
Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850


DiamoDent®

Compatible with: - 6.0mm (5.7mm Platform) Screw-Vent Zimmer® implants
- 5.7mm Diameter (5.7mm Platform) Legacy® Implant Direct® implants

UCLA ABUTMENT, Gold/Clear Plastic, Yellow

S110 Hexed, 0.85mm Collar, 5.7mm Profile, TiteFit™
S101 Non-Hexed, 0.85mm Collar, 5.7mm Profile

Requires ComScrew catalog number P655



S110 S101

CEMENT-ON POST, Titanium

S505 1mm Collar, 6.5mm Profile, TiteFit™
S510 2mm Collar, 6.5mm Profile, TiteFit™
S515 4mm Collar, 6.5mm Profile, TiteFit™
S535 3mm Collar, 6.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



S505 S510 S515 S535

ComScrew™, Titanium

P655 .050" Hex
30 Ncm



GUIDE PIN/WAXING SCREW, Titanium

P680 .050" Hex



THREADED IMPRESSION COPING, Titanium, Yellow

S700 Closed Tray, 6mm Profile



DíamoDent®

Compatible with: - 6.0mm (5.7mm Platform) Screw-Vent Zimmer® implants
- 5.7mm Diameter (5.7mm Platform) Legacy® Implant Direct® implants

THREADED IMPRESSION COPING, Titanium, Yellow

S750

Open Tray, 6.4mm Profile



TEMPORARY ABUTMENT, Titanium

S800
S801

Hexed, 1mm Collar, 6.4mm Profile
Non-Hexed, 1mm Collar, 6.4mm Profile



Requires ComScrew catalog number P655

S800

S801

IMPLANT REPLICA, Aluminum, Yellow

S910




ZIRCONIA ABUTMENT BASE, Titanium, Yellow

S1000 Hexed, 0.35 Collar, 6.2mm Profile
Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850


DiamoDent®

Compatible with: 3.5mm (3.5mm Platform) BioHorizon® implants

UCLA ABUTMENT, Gold/Clear Plastic, Yellow

AP110 Hexed, 0.8mm Collar, 4.5mm, Profile, TiteFit™
AP101 Non-Hexed, 0.8mm Collar, 4.5mm Profile

Requires ComScrew catalog number P655



AP110

AP101

UCLA ABUTMENT, Clear Plastic, Yellow

AP200 Hexed, 0.9mm Collar, 4.3mm Profile
AP201 Non-Hexed, 0.9mm Collar, 4.3mm Profile

Requires ComScrew catalog number P655



AP200

AP201

OVERDENTURE ABUTMENT, Clear Plastic, Yellow

AP401 Non-Hexed

Requires ComScrew catalog number P655

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

AP505 1mm Collar, 4.5mm Profile, TiteFit™
AP510 2mm Collar, 4.5mm Profile, TiteFit™
AP515 4mm Collar, 4.5mm Profile, TiteFit™
AP520 2mm Collar, 5.5mm Profile, TiteFit™
AP525 4mm Collar, 5.5mm Profile, TiteFit™
AP535 3mm Collar, 5.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



AP505

AP510

AP515

AP520

AP525

AP535

ComScrew®, Titanium

P655 .050" Hex
30 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 3.5mm (3.5mm Platform) BioHorizon® implants

GUIDE PIN/WAXING SCREW, Titanium

P680

.050" Hex



THREADED IMPRESSION COPING, Titanium, Yellow

AP700

Closed Tray, 4.6mm Profile



THREADED IMPRESSION COPING, Titanium, Yellow

AP750

Open Tray, 4.7mm Profile



TEMPORARY ABUTMENT, Titanium

AP800

Hexed, 1mm Collar, 4.5mm Profile

AP801

Non-Hexed, 1mm Collar, 4.5mm Profile



Requires ComScrew catalog number P655

AP800

AP801

MPLANT REPLICA, Aluminum, Yellow

AP910



ZIRCONIA ABUTMENT BASE, Titanium, Yellow

AP1000


Hexed, 0.35 Collar, 4.0mm Profile

Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850

DiamoDent®

Compatible with: 4.0mm (4.5mm Platform) BioHorizon® implants

UCLA ABUTMENT, Gold/Clear Plastic, Green

AR110 Hexed, 0.8mm Collar, 4.7mm Profile, TiteFit™
AR101 Non-Hexed, 0.8mm Collar, 4.7mm Profile

Requires ComScrew catalog number P655



AR110



AR101

UCLA ABUTMENT, Clear Plastic, Green

AR200 Hexed, 0.7mm Collar, 4.6mm Profile
AR201 Non-Hexed, 0.7mm Collar, 4.6mm Profile

Requires ComScrew catalog number P655



AR200



AR201

OVERDENTURE ABUTMENT, Clear Plastic, Green

AR401 Non-Hexed

Requires ComScrew catalog number P655

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

AR505 1mm Collar, 5.5mm Profile, TiteFit™
AR510 2mm Collar, 5.5mm Profile, TiteFit™
AR515 4mm Collar, 5.5mm Profile, TiteFit™
AR520 2mm Collar, 6.5mm Profile, TiteFit™
AR525 4mm Collar, 6.5mm Profile, TiteFit™
AR535 3mm Collar, 6.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



AR505



AR510



AR515



AR520



AR525



AR535

ComScrew™, Titanium

P655 .050" Hex
30 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.0mm (4.5mm Platform) BioHorizon® implants

GUIDE PIN/WAXING SCREW, Titanium

P680

.050" Hex



THREADED IMPRESSION COPING, Titanium, Green

AR700

Closed Tray, 4.9mm Profile



THREADED IMPRESSION COPING, Titanium, Green

AR750

Open Tray, 5.2mm Profile



TEMPORARY ABUTMENT, Titanium

AR800

Hexed, 1mm Collar, 5.4mm Profile

AR801

Non-Hexed, 1mm Collar, 5.4mm Profile



AR800

AR801

Requires ComScrew catalog number P655

IMPLANT REPLICA, Aluminum, Green

AR910



ZIRCONIA ABUTMENT BASE, Titanium, Green

AR1000

Hexed, 0.35 Collar, 5.0mm Profile

Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 5.0mm & 6.0mm (5.7mm Platform) BioHorizon® implants

UCLA ABUTMENT, Gold/Clear Plastic, Blue

AS110	Hexed, 0.85mm Collar, 5.7mm Profile, TiteFit™
AS101	Non-Hexed, 0.85mm Collar, 5.7mm Profile

Requires ComScrew catalog number P655



AS110



AS101

CEMENT-ON POST, Titanium

AS505	1mm Collar, 6.5mm Profile, TiteFit™
AS510	2mm Collar, 6.5mm Profile, TiteFit™
AS515	4mm Collar, 6.5mm Profile, TiteFit™
AS535	3mm Collar, 6.5mm Profile, TiteFit™

Requires ComScrew catalog number P655



AS505



AS510



AS515



AS535

ComScrew™, Titanium

P655 .050" Hex

30 Ncm



GUIDE PIN/WAXING SCREW, Titanium

P680 .050" Hex



THREADED IMPRESSION COPING, Titanium, Blue

AS700 Closed Tray, 6mm Profile



DíamoDent®

Compatible with: 5.0mm & 6.0mm (5.7mm Platform) BioHorizon® implants

THREADED IMPRESSION COPING, Titanium, Blue

AS750 Open Tray, 6.4mm Profile



TEMPORARY ABUTMENT, Titanium

AS800 Hexed, 1mm Collar, 6.4mm Profile
AS801 Non-Hexed, 1mm Collar, 6.4mm Profile

Requires ComScrew catalog number P655



AS800



AS801

IMPLANT REPLICA, Aluminum, Blue

AS910




ZIRCONIA ABUTMENT BASE, Titanium, Blue

AS1000 Hexed, 0.35 Collar, 6.2mm Profile
Requires ComScrew catalog number P655



SCREW AND ABUTMENT REMOVAL TOOL

T5190



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: : 3.4mm Platform Internal BIOMET® 3i® Certain implants

UCLA ABUTMENT, Gold/Clear Plastic, Purple

AJ100 Hexed, 0.50mm Collar, 3.6mm Profile
Requires Small Diameter ComScrew™

AJ101 Non-Hexed, 0.50mm Collar, 3.6mm Profile
Requires Large Diameter ComScrew™



AJ100



AJ101

UCLA ABUTMENT, Clear Plastic, Purple

AJ201 Non-Hexed, 1.0mm Collar, 4.2mm Profile
Requires Large Diameter ComScrew™



AJ201

CEMENT-ON POST, Titanium

AJ505 1mm Collar, 4.4mm Profile
AJ510 2mm Collar, 4.4mm Profile
AJ535 3mm Collar, 4.4mm Profile

Requires Small Diameter ComScrew™



AJ505



AJ510



AJ535

SMALL DIAMETER ComScrew™

AJ600 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ650 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ600



AJ650

LARGE DIAMETER ComScrew™

AJ601 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ651 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ601



AJ651

CUSTOMER SERVICE
888-281-8850

DiamoDent®

Compatible with: 3.4mm Platform Internal BIOMET® 3i® Certain implants

GUIDE PIN/WAXING SCREW, Titanium

AJ680 Small Diameter, 1.2mm (.048") Hex
AJ681 Large Diameter, 1.2mm (.048") Hex



THREADED IMPRESSION COPING, Titanium, Purple

AJ700 Closed Tray, 4.0mm Profile



THREADED IMPRESSION COPING, Titanium, Purple

AJ750 Open Tray, 4.1mm Profile



TEMPORARY ABUTMENT, Titanium

AJ800 Hexed, 1.0mm Collar, 4.2mm Profile
Requires Small Diameter ComScrew™

AJ801 Non-Hexed, 1.0mm Collar, 4.2mm Profile
Requires Large Diameter ComScrew™



AJ800



AJ801

IMPLANT REPLICA, Aluminum, Purple

AJ910



ZIRCONIA ABUTMENT BASE, Titanium, Purple

AJ1000 Hexed, 0.35 Collar, 4.0mm Profile
Requires Small Diameter ComScrew™



DíamoDent®

Compatible with: : 4.1mm Platform Internal BIOMET® 3i® Certain implants

UCLA ABUTMENT, Gold/Clear Plastic, Blue

AK100 Hexed, 0.50mm Collar, 4.5mm Profile
Requires Small Diameter ComScrew™

AK101 Non-Hexed, 0.50mm Collar, 4.5mm Profile
Requires Large Diameter ComScrew™



AK100



AK101

UCLA ABUTMENT, Clear Plastic, Blue

AK200 Hexed, 1.0mm Collar, 4.8mm Profile
Requires Small Diameter ComScrew™

AK201 Non-Hexed, 1.0mm Collar, 4.8mm Profile
Requires Large Diameter ComScrew™



AK200



AK201

CEMENT-ON POST, Titanium

AK505 1mm Collar, 5.1mm Profile
AK510 2mm Collar, 5.1mm Profile
AK515 4mm Collar, 5.1mm Profile
AK520 2mm Collar, 5.9mm Profile
AK525 4mm Collar, 5.9mm Profile

Requires Small Diameter ComScrew™



AK505



AK510



AK515



AK520



AK525

SMALL DIAMETER ComScrew™

AJ600 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ650 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ600



AJ650

LARGE DIAMETER ComScrew™

AJ601 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ651 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ601



AJ651

CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.1mm Platform Internal BIOMET® 3i® Certain implants

Guide Pin/Waxing Screw, Titanium

AJ680 Small Diameter, 1.2mm (.048") Hex
AJ681 Large Diameter, 1.2mm (.048") Hex



THREADED IMPRESSION COPING, Titanium, Blue

AK700 Closed Tray, 4.3mm Profile



THREADED IMPRESSION COPING, Titanium, Blue

AK750 Open Tray, 4.5mm Profile



TEMPORARY ABUTMENT, Titanium

AK800 Hexed, 1.0mm Collar, 4.9mm Profile
Requires Small Diameter ComScrew™

AK801 Non-Hexed, 1.0mm Collar, 4.9mm Profile
Requires Large Diameter ComScrew™



AK800



AK801

IMPLANT REPLICA, Aluminum, Blue

AK910



ZIRCONIA ABUTMENT BASE, Titanium, Blue

AK1000 Hexed, 0.35 Collar, 4.6mm Profile
Requires Small Diameter ComScrew™



DiamoDent®

Compatible with: : 5.0mm Platform Internal BIOMET® 3i® Certain implants

UCLA ABUTMENT, Gold/Clear Plastic, Yellow

AL100 Hexed, 0.50mm Collar, 5.5mm Profile
Requires Small Diameter ComScrew™

AL101 Non-Hexed, 0.50mm Collar, 5.5mm Profile
Requires Large Diameter ComScrew™



AL100



AL101

UCLA ABUTMENT, Clear Plastic, Yellow

AL200 Hexed, 1.0mm Collar, 5.5mm Profile
Requires Small Diameter ComScrew™

AL201 Non-Hexed, 1.0mm Collar, 5.5mm Profile
Requires Large Diameter ComScrew™



AL200



AL201

CEMENT-ON POST, Titanium

AL505 1mm Collar, 6.0mm Profile
AL510 2mm Collar, 6.0mm Profile
AL515 4mm Collar, 6.0mm Profile
AL520 2mm Collar, 6.80mm Profile
AL525 4mm Collar, 6.80mm Profile

Requires Small Diameter ComScrew™



AL505



AL510



AL515



AL520



AL525

SMALL DIAMETER ComScrew™

AJ600 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ650 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ600



AJ650

LARGE DIAMETER ComScrew™

AJ601 Gold, 1.2mm (.048") Hex, 20Ncm
Gold alloy

AJ651 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ601



AJ651

CUSTOMER SERVICE
888-281-8850

DiamoDent®

Compatible with: 5.0mm Platform Internal BIOMET® 3i® Certain implants

GUIDE PIN/WAXING SCREW, Titanium

AJ680 Small Diameter, 1.2mm (.048") Hex
AJ681 Large Diameter, 1.2mm (.048") Hex



THREADED IMPRESSION COPING, Titanium, Yellow

AL700 Closed Tray, 5.5mm Profile



THREADED IMPRESSION COPING, Titanium, Yellow

AL750 Open Tray, 5.5mm Profile



TEMPORARY ABUTMENT, Titanium

AL800 Hexed, 1.0mm Collar, 5.8mm Profile
Requires Small Diameter ComScrew™

AL801 Non-Hexed, 1.0mm Collar, 5.8mm Profile
Requires Large Diameter ComScrew™



AL800

AL801

IMPLANT REPLICA, Aluminum, Yellow

AL910



ZIRCONIA ABUTMENT BASE, Titanium, Yellow

AL1000 Hexed, 0.35 Collar, 5.5mm Profile
Requires Small Diameter ComScrew™



DíamoDent®

Compatible with: : 6.0mm Platform Internal BIOMET® 3i® Certain implants

UCLA ABUTMENT, Gold/Clear Plastic, Green

AM100 Hexed, 0.50mm Collar, 6.5mm Profile
Requires Small Diameter ComScrew™

AM101 Non-Hexed, 0.50mm Collar, 6.5mm Profile
Requires Large Diameter ComScrew™



AM100



AM101

UCLA ABUTMENT, Clear Plastic, Green

AM200 Hexed, 1.0mm Collar, 6.5mm Profile
Requires Small Diameter ComScrew™

AM201 Non-Hexed, 1.0mm Collar, 6.5mm Profile
Requires Large Diameter ComScrew™



AM200



AM201

CEMENT-ON POST, Titanium

AM505 1mm Collar, 7.0mm Profile
AM510 2mm Collar, 7.0mm Profile
AM515 4mm Collar, 7.0mm Profile

Requires Small Diameter ComScrew™



AM505



AM510



AM515

SMALL DIAMETER ComScrew™

AJ600 Gold alloy Gold, 1.2mm (.048") Hex, 20Ncm

AJ650 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ600



AJ650

CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 6.0mm Platform Internal BIOMET® 3i® Certain implants

LARGE DIAMETER ComScrew™

AJ601 Gold, 1.2mm (.048") Hex, 20Ncm
AJ651 Titanium, 1.2mm (.048") Hex, 20Ncm



AJ601



AJ651

GUIDE PIN/WAXING SCREW, Titanium

AJ680 Small Diameter, 1.2mm (.048") Hex
AJ681 Large Diameter, 1.2mm (.048") Hex



AJ680



AJ681

THREADED IMPRESSION COPING, Titanium, Green

AM700 Closed Tray, 6.5mm Profile



IMPLANT REPLICA, Aluminum, Green

AM910



ZIRCONIA ABUTMENT BASE, Titanium, Green

AM1000 Hexed, 0.35 Collar, 6.5mm Profile
Requires Small Diameter ComScrew™



DiamoDent®

Compatible with: 4.1mm / 4.8mm Regular Neck (RN) ITI Straumann® implants

UCLA ABUTMENT, Gold/Clear Plastic, Red

Y100 Engaging, 0.4mm Collar, 5mm Profile
Y101 Non-Engaging, 0.4mm Collar, 5mm Profile

Requires ComScrew catalog number Y665



Y100



Y101

SOLID ABUTMENT PLASTIC COPING

Y300 Engaging
Y301 Non-Engaging



UCLA COPING, Gold/Clear Plastic, Red

Y400 Engaging, 0.2mm Collar, 5mm Profile
Y401 Non-Engaging, 0.2mm Collar, 5mm Profile

Requires ComScrew catalog number Y670



CEMENT-ON POST, Titanium

Y500 1mm Collar, 5.75mm Profile

Requires ComScrew catalog number Y665



SOLID ABUTMENT®, Titanium

Y550 4mm Height, Yellow
Y555 5.5mm Height, Gray
Y560 7mm Height, Blue



Y550



Y555



Y560

ABUTMENT SCREW, Titanium

Y665 SCS®
35 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.1mm / 4.8mm Regular Neck (RN) ITI Straumann® implants

COPING SCREW, Titanium

Y670 SCS®
35 Ncm



THREADED IMPRESSION COPING, Red

Y700 Closed Tray, 5.4mm Profile



THREADED IMPRESSION COPING, Red

Y750 Open Tray, 5.4mm Profile



TEMPORARY ABUTMENT, Titanium

Y800 Engaging, 0.3mm Collar, 5mm Profile
Y801 Non-Engaging, 0.3mm Collar, 5mm Profile



Requires ComScrew catalog number Y665

IMPLANT REPLICA, Stainless Steel

Y900



SOLID ABUTMENT® ANALOG, Aluminum

Y950 4mm Height, Yellow
Y955 5.5mm Height, Gray
Y960 7mm Height, Blue



Y950



Y955



Y960


CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.8mm Wide Neck (WN) ITI Straumann® implants

UCLA ABUTMENT, Gold/Clear Plastic

Z100 Engaging, 0.4mm Collar, 6.7mm Profile
Z101 Non-Engaging, 0.4mm Collar, 6.7mm Profile



Z100

Z101

Requires ComScrew catalog number Y665

SOLID ABUTMENT® PLASTIC COPING

Z300 Engaging
Z301 Non-Engaging



CEMENT-ON POST, Titanium

Z500 1.25mm Collar, 7.5mm Profile

Requires ComScrew catalog number Y665



SOLID ABUTMENT®, Titanium

Z550 4mm Height, Green, SCS
Z555 5.5mm Height, Brown, SCS



Z550



Z555

ABUTMENT SCREW, Titanium

Y665 SCS®

35 Ncm



THREADED IMPRESSION COPING, Titanium

Z700 Closed Tray, 7mm Profile



DíamoDent®

Compatible with: 4.8mm Wide Neck (WN) ITI Straumann® implants

THREADED IMPRESSION COPING, Titanium

Z750 Open Tray, 7mm Profile



TEMPORARY ABUTMENT, Titanium

Z800 Engaging, 0.75mm Collar, 7mm Profile
Z801 Non-Engaging, 0.75mm Collar, 7mm Profile

Requires ComScrew catalog number Y665



IMPLANT REPLICA, Stainless Steel

Z900

Suitable for soldering



SOLID ABUTMENT® ANALOG, Aluminum

Z950 4mm Height, Green
Z955 5.5mm Height, Brown



Z950



Z955

DíamoDent®

Compatible with: 3.3mm Narrow Neck (NN) ITI Straumann® implants

UCLA ABUTMENT, Gold/Clear Plastic

X100	Engaging, 0.4mm Collar, 3.8mm Profile
X101	Non-Engaging, 0.4mm Collar, 3.8mm Profile

Requires ComScrew catalog number X665



ComScrew™, Titanium, Yellow

X665	SCS®
------	------

35 Ncm



THREADED IMPRESSION COPING, Titanium

X750	Open Tray, 4mm Profile
------	------------------------



TEMPORARY ABUTMENT, Titanium

X800	Engaging, 1mm Collar, 4mm Profile
X801	Non-Engaging, 1mm Collar, 4mm Profile

Requires ComScrew catalog number X665



IMPLANT REPLICA, Stainless Steel

X900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 3.4mm Frialit-2 & XiVE Friadent implants®

UCLA ABUTMENT, Gold/Clear Plastic

AE100 Hexed, 0.8mm Collar, 3.6mm Profile

Requires ComScrew catalog number AE650



CEMENT-ON POST, Titanium

AE500 2mm Collar, 4.2mm Profile
AE510 4mm Collar, 4.2mm Profile

Requires ComScrew catalog number AE650



AE500



AE510

ComScrew™, Titanium

AE650 1.2mm (.048") Hex

20 Ncm



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 3.8mm Frialit-2 & XiVE Friadent implants®

UCLA ABUTMENT, Gold/Clear Plastic, Yellow

AF100 Hexed, 0.8mm Collar, 4mm Profile

Requires ComScrew catalog number AE650



CEMENT-ON POST, Titanium

AF500 2mm Collar, 4.5mm Profile
AF510 4mm Collar, 4.5mm Profile

Requires ComScrew catalog number AE650



AF500



AF510

ComScrew™, Titanium

AE650 1.2mm (.048") Hex

20 Ncm



IMPLANT REPLICA, Stainless Steel

AF900

Suitable for soldering



IMPLANT REPLICA, Aluminum, Yellow

AF910



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.5mm Frialit-2 & XiVE Friadent implants®

CEMENT-ON POST, Titanium

AG500 2mm Collar, 5mm Profile
AG510 4mm Collar, 5mm Profile

Requires ComScrew catalog number AE650



AG500



AG510

ComScrew™, Titanium

AE650 1.2mm (.048") Hex

20 Ncm



IMPLANT REPLICA, Aluminum, Blue

AG910



DíamoDent®

Compatible with: 5.5mm Frialit-2 & XiVE Friadent implants®

UCLA ABUTMENT, Gold/Clear Plastic, Red

AH100 Hexed, 0.8mm Collar, 5.6mm Profile

Requires ComScrew catalog number AE650



CEMENT-ON POST, Titanium

AH510 4mm Collar, 6mm Profile

Requires ComScrew catalog number AE650



AH510

ComScrew™, Titanium

AE650 1.2mm (.048") Hex

20 Ncm



IMPLANT REPLICA, Stainless Steel

AH900

Suitable for soldering



IMPLANT REPLICA, Aluminum, Red

AH910



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 3.25mm & 3.3mm MicroMiniplant BIOMET® 3i® implants

UCLA ABUTMENT, Gold/Clear Plastic

J100 Hexed, 0.85mm Collar, 3.6mm Profile
J101 Non-Hexed, 0.85mm Collar, 3.6mm Profile

Requires ComScrew



CEMENT-ON POST, Titanium

J500 2mm Collar, 4mm Profile
J510 4mm Collar, 4mm Profile

Requires ComScrew



J500



J510

ComScrew™, Gold

A600 Square

32 Ncm

Gold alloy



ComScrew™, Titanium

A650 1.2mm (.048") Hex
A655 .050" Hex

20 Ncm



A650



A655

THREADED IMPRESSION COPING, Titanium, Purple

J750 Open Tray, 3.8mm Profile



IMPLANT REPLICA, Stainless Steel

J900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: -3.75 & 4mm Regular Platform (RP) Branemark® implants
-Osseotite® XP 3/4, 3.25 & 3.3mm Miniplants, 3.75 & 4mm Standard
BIOMET® 3i® implants
-3.8 & 4.5mm HL Steri-Oss® implants
-3.75, 4.0 & 4.2mm Regular Diameter (RD) RESTORE® Lifecore®
/ Keystone® implants

UCLA ABUTMENT, Gold/Clear Plastic

A100	Hexed, 0.8mm Collar, 4.5mm Profile
A101	Non-Hexed, 0.8mm Collar, 4.5mm Profile

Requires ComScrew



UCLA ABUTMENT, Clear Plastic

A200	Hexed, 1.2mm Collar, 4.4mm Profile
A201	Non-Hexed, 1.2mm Collar, 4.4mm Profile

Requires ComScrew



OVERDENTURE ABUTMENT, Clear Plastic

A401	Non-Hexed
------	-----------

Requires ComScrew



For use with screw retained dentures, Hader Bar and distal extension attachment cases.

OVERDENTURE ABUTMENT, Titanium

A451	Non-Hexed
------	-----------

Requires ComScrew

Suitable for laser welding to Titanium Hader Bar



CEMENT-ON POST, Titanium

A500	2mm Collar, 5mm Profile
A510	4mm Collar, 5mm Profile

Requires ComScrew



A500



A510

CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with:

- 3.75 & 4mm Regular Platform (RP) Branemark® implants
- Osseotite® XP 3/4, 3.25 & 3.3mm Miniplants, 3.75 & 4mm Standard BIOMET® 3i® implants
- 3.8 & 4.5mm HL Steri-Oss® implants
- 3.75, 4.0 & 4.2mm Regular Diameter (RD) RESTORE® Lifecore® / Keystone® implants

ComScrew™, Gold

A600 Square

32Ncm

Gold alloy



Try-In ComScrew™, Titanium, Blue

A640 Square



ComScrew™, Titanium

A650 1.2mm (.048")Hex, 20 Ncm
A655 .050" Hex, 20 Ncm
A660 UniGrip®, 35 Ncm



A650



A655



A660

THREADED IMPRESSION COPING, Titanium

A700 Closed Tray, 4.6mm Profile



TEMPORARY ABUTMENT, Titanium

A800 Hexed, 1.4mm Collar, 4.5mm Profile
A801 Non-Hexed, 1.4mm Collar, 4.5mm Profile



Requires ComScrew

IMPLANT REPLICA, Stainless Steel

A900

Suitable for soldering



DíamoDent®

Compatible with: 5.0mm BIOMET® 3i® implants

UCLA ABUTMENT, Gold/Clear Plastic

C100 Hexed, 0.8mm Collar, 5.5mm Profile
C101 Non-Hexed, 0.8mm Collar, 5.5mm Profile

Requires ComScrew



UCLA ABUTMENT, Clear Plastic

C200 Hexed, 1.3mm Collar, 5.4mm Profile
C201 Non-Hexed, 1.3mm Collar, 5.4mm Profile

Requires ComScrew



OVERDENTURE ABUTMENT, Clear Plastic

C401 Non-Hexed

Requires ComScrew

For use with screw retained dentures, Hader Bar and distal extension attachment cases.



CEMENT-ON POST, Titanium

C500 2mm Collar, 6mm Profile
C510 4mm Collar, 6mm Profile

Requires ComScrew



C500



C510

ComScrew™, Gold

A600 Square

32 Ncm

Gold alloy



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 5.0mm BIOMET® 3i® implants

Try-In ComScrew™, Titanium , Blue

A640

Square



ComScrew™, Titanium

A650

1.2mm (.048") Hex

A655

.050" Hex

20 Ncm



A650



A655

THREADED IMPRESSION COPING, Titanium

C700

Closed Tray, 6mm Profile



THREADED IMPRESSION COPING, Titanium

C750

Open Tray, 6mm Profile



IMPLANT REPLICA, Stainless Steel

C900

Suitable for soldering



DíamoDent®

Compatible with: 6.0mm BIOMET® 3i® implants

UCLA ABUTMENT, Gold/Clear Plastic

D100 Hexed, 0.85mm Collar, 6.5mm Profile

Requires ComScrew



CEMENT-ON POST, Titanium

D500 2mm Collar, 7mm Profile

D510 4mm Collar, 7mm Profile

Requires ComScrew



D500



D510

ComScrew™, Gold

A600 Square

32 Ncm
Gold alloy



Try-In ComScrew™, Titanium , Blue

A640 Square



ComScrew™, Titanium

A650 1.2mm (.048") Hex

A655 .050" Hex

20 Ncm



A650



A655

THREADED IMPRESSION COPING, Titanium

D700 Closed Tray, 7mm Profile



THREADED IMPRESSION COPING, Titanium

D750 Open Tray, 7mm Profile



IMPLANT REPLICA, Stainless Steel

D900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

Compatible with: 5.0mm Wide Platform (WP) Branemark® implants

UCLA ABUTMENT, Gold/Clear Plastic

B100 Hexed, 0.85mm Collar, 5.5mm Profile
B101 Non-Hexed, 0.85mm Collar, 5.5mm Profile

Requires ComScrew



UCLA ABUTMENT, Clear Plastic

B201 Non-Hexed, 1.3mm Collar, 5.4mm Profile

Requires ComScrew



OVERDENTURE ABUTMENT, Clear Plastic

B401 Non-Hexed

Requires ComScrew



For use with screw retained dentures, Hader Bar and distal extension attachment cases.

CEMENT-ON POST, Titanium

B500 1.4mm Collar, 6mm Profile

Requires ComScrew



ComScrew™, Gold

B600 Square

32 Ncm
Gold alloy



ComScrew™, Titanium

N650 1.2mm (.048") Hex, 20 Ncm
B660 UniGrip®, 35 Ncm



N650



B660

IMPLANT REPLICA, Stainless Steel

B900

Suitable for soldering



DíamoDent®

Compatible with: **3.3mm & 3.4mm Small Diameter (SD) RESTORE® Lifecore® / Keystone® implants**

UCLA ABUTMENT, Gold/Clear Plastic

E100	Hexed, 0.8mm Collar, 3.6mm Profile
E101	Non-Hexed, 0.8mm Collar, 3.6mm Profile

Requires ComScrew



CEMENT-ON POST, Titanium

E500	2mm Collar, 4mm Profile
E510	4mm Collar, 4mm Profile

Requires ComScrew



E500



E510

ComScrew™, Titanium

E650	1.2mm (.048") Hex
E655	.050" Hex

20 Ncm



E650



E655

THREADED IMPRESSION COPING, Titanium

E750	Open Tray, 4mm Profile
------	------------------------



IMPLANT REPLICA, Stainless Steel

E900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

DiamoDent®

Compatible with: : 5.0mm & 6.0mm Wide Diameter (WD) RESTORE®
Lifecore® / Keystone® implants

UCLA ABUTMENT, Gold/Clear Plastic

N101 Non-Hexed, 0.85mm Collar, 5.5mm Profile

Requires ComScrew



CEMENT-ON POST, Titanium

N500 2mm Collar, 6mm Profile
N510 4mm Collar, 6mm Profile

Requires ComScrew



N500



N510

ComScrew™, Gold

B600 Square

32 Ncm
Gold alloy



ComScrew™, Titanium

N650 1.2mm (.048") Hex

20 Ncm



N650

THREADED IMPRESSION COPING, Titanium

N700 Closed Tray, 5.5mm Profile



THREADED IMPRESSION COPING, Titanium

N750 Open Tray, 5.8mm Profile



IMPLANT REPLICA, Stainless Steel

N900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 4.3mm External Hex Replace Steri-Oss® implants

UCLA ABUTMENT, Clear Plastic

H200	Hexed, 1.5mm Collar, 4.6mm Profile
H201	Non-Hexed, 1.5mm Collar, 4.6mm Profile

Requires ComScrew



CEMENT-ON POST, Titanium

H500	1.4mm Collar, 5.2mm Profile
------	-----------------------------

Requires ComScrew



ComScrew™, Gold

A600	Square
------	--------

32 Ncm
Gold alloy



ComScrew™, Titanium

A650	1.2mm (.048") Hex, 20 Ncm
A655	.050" Hex, 20 Ncm
A660	UniGrip®, 35 Ncm



A650



A655



A660

IMPLANT REPLICA, Stainless Steel

H900

Suitable for soldering



CUSTOMER SERVICE
888-281-8850

DíamoDent®

Compatible with: 5.0mm External Hex Replace® and HL Threaded Steri-Oss® implants

UCLA ABUTMENT, Gold/Clear Plastic

K101 Non-Hexed, 0.8mm Collar, 5.5mm Profile

Requires ComScrew



UCLA ABUTMENT, Clear Plastic

K200 Hexed, 1.3mm Collar, 5.4mm Profile
K201 Non-Hexed, 1.3mm Collar, 5.4mm Profile

Requires ComScrew



OVERDENTURE ABUTMENT, Clear Plastic

K401 Non-Hexed

Requires ComScrew



For use with screw retained dentures, Hader Bar and distal extension attachment cases.

CEMENT-ON POST, Titanium

K500 1.4mm Collar, 6mm Profile

Requires ComScrew



ComScrew™, Gold

A600 Square

32 Ncm
Gold alloy



ComScrew™, Titanium

A650 1.2mm (.048") Hex, 20 Ncm
A655 .050" Hex, 20 Ncm
A660 UniGrip®, 35 Ncm



A650



A655



A660

IMPLANT REPLICA, Stainless Steel

K900





Suitable for soldering



CUSTOMER SERVICE
888-281-8850

Tools and Accessories

SCREWDRIVERS

T5100	.035" Hex, 28.2mm Length	
T5105	.035" Hex, 21mm Length	
T5120	1.2mm (.048") Hex, 28.2mm Length	
T5125	1.2mm (.048") Hex, 21mm Length	
T5140	.050" Hex, 28.2mm Length	
T5145	.050" Hex, 21mm Length	
T5160	Square, 28.2mm Length	
T5170	UniGrip®, 28.2mm Length	
T5175	UniGrip®, 21mm Length	
T5180	SCS®, 28.2mm Length	
T5185	SCS®, 21mm Length	

SCREW AND ABUTMENT REMOVAL TOOL*

T5190*
*U.S. and International patents pending



SCREW, Gold

T5220 Slot




SCREW, Titanium

T5240 Slot



HADER BAR, Clear Plastic

T5320-6 3mm H x 50mm L, 6/Pkg
T5340-6 8mm H x 50mm L, 6/Pkg




Compatible with all retention clips


Tools and Accessories

DRIVER TIPS


T5400 1.2mm (.048) Hex, Long




T5405 1.2mm (.048) Hex, Short




T5410 .050" Hex, Long




T5415 .050" Hex, Short




T5420 Square, Long




T5425 Square, Short




T5430 Unigrip®, Long




T5435 Unigrip®, Short



T5440 SCS®, Long



T5445 SCS®, Short



DíamoDent®

Tools and Accessories

LAB SCREWDRIVER

T5500	1.2mm (.048") Hex	
T5510	.050" Hex	
T5520	Square	
T5530	UniGrip®	
T5540	SCS®	
T5550	Slotted	

LAB SCREWDRIVER HOLDER

T5580	Holder
-------	--------

LAB SCREWDRIVER KIT

T5590

Contains six Lab Screwdrivers,
Holder and two mounting screws



CUSTOMER SERVICE
888-281-8850



Price List

Effective March 2013

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
A100	\$90	\$88.50	\$87	\$85	\$85
A101	\$90	\$88.50	\$87	\$85	\$85
A200	\$10	\$9.75	\$9.50	\$9.25	\$9
A201	\$10	\$9.75	\$9.50	\$9.25	\$9
A401	\$10	\$9.75	\$9.50	\$9.25	\$9
A451	\$25	\$24.50	\$24	\$23.50	\$23
A500	\$50	\$48	\$46.50	\$46	\$45.50
A510	\$50	\$48	\$46.50	\$46	\$45.50
A600	\$50	\$49	\$48	\$47	\$47
A640	\$12	\$11.75	\$11.50	\$11	\$10.50
A650	\$12	\$11.75	\$11.50	\$11	\$10.50
A655	\$12	\$11.75	\$11.50	\$11	\$10.50
A660	\$14	\$13.50	\$13	\$12.50	\$12
A700	\$35	\$34	\$33	\$33	\$33
A800	\$25	\$24.50	\$24	\$23.50	\$23
A801	\$25	\$24.50	\$24	\$23.50	\$23
A900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
B100	\$100	\$99	\$97	\$95	\$95
B101	\$100	\$99	\$97	\$95	\$95
B201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
B401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
B500	\$50	\$48	\$46.50	\$46	\$45.50
B600	\$55	\$54	\$53	\$52	\$51
B660	\$14	\$13.50	\$13	\$12.50	\$12
B900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
C100	\$100	\$99	\$97	\$95	\$95
C101	\$100	\$99	\$97	\$95	\$95
C200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
C201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
C401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
C500	\$50	\$48	\$46.50	\$46	\$45.50
C510	\$50	\$48	\$46.50	\$46	\$45.50
C700	\$35	\$34	\$33	\$33	\$33
C750	\$35	\$34	\$33	\$33	\$33
C900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
D100	\$110	\$108	\$106	\$104	\$104
D500	\$50	\$48	\$46.50	\$46	\$45.50
D510	\$50	\$48	\$46.50	\$46	\$45.50
D700	\$35	\$34	\$33	\$33	\$33
D750	\$35	\$34	\$33	\$33	\$33
D900	\$12	\$11.75	\$11.50	\$11	\$10.50
E100	\$90	\$88.50	\$87	\$85	\$85
E101	\$90	\$88.50	\$87	\$85	\$85
E500	\$50	\$48	\$46.50	\$46	\$45.50
E510	\$50	\$48	\$46.50	\$46	\$45.50
E650	\$12	\$11.75	\$11.50	\$11	\$10.50
E655	\$12	\$11.75	\$11.50	\$11	\$10.50
E750	\$35	\$34	\$33	\$33	\$33
E900	\$10	\$9.60	\$9.20	\$8.85	\$8.50

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
H200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
H201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
H500	\$50	\$48	\$46.50	\$46	\$45.50
H900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
J100	\$90	\$88.50	\$87	\$85	\$85
J101	\$90	\$88.50	\$87	\$85	\$85
J500	\$50	\$48	\$46.50	\$46	\$45.50
J510	\$50	\$48	\$46.50	\$46	\$45.50
J750	\$35	\$34	\$33	\$33	\$33
J900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
K101	\$100	\$99	\$97	\$95	\$95
K200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
K201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
K401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
K500	\$50	\$48	\$46.50	\$46	\$45.50
K900	\$10	\$9.60	\$9.20	\$8.85	\$8.50
N101	\$100	\$99	\$97	\$95	\$95
N500	\$50	\$48	\$46.50	\$46	\$45.50
N510	\$50	\$48	\$46.50	\$46	\$45.50
N650	\$12	\$11.75	\$11.50	\$11	\$10.50
N700	\$35	\$34	\$33	\$33	\$33
N750	\$35	\$34	\$33	\$33	\$33
N900	\$12	\$11.75	\$11.50	\$11	\$10.50
P101	\$80	\$78	\$76	\$74	\$74
P110	\$80	\$78	\$76	\$74	\$74
P200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
P201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
P401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
P505	\$55	\$53	\$51	\$50	\$49
P510	\$55	\$53	\$51	\$50	\$49
P515	\$55	\$53	\$51	\$50	\$49
P520	\$55	\$53	\$51	\$50	\$49
P525	\$55	\$53	\$51	\$50	\$49
P535	\$55	\$53	\$51	\$50	\$49
P655	\$12	\$11.75	\$11.50	\$11	\$10.50
P680	\$13	\$12.75	\$12.50	\$12	\$11.50
P700	\$35	\$34	\$33	\$33	\$33
P750	\$35	\$34	\$33	\$33	\$33
P800	\$30	\$29	\$28	\$28	\$28
P801	\$30	\$29	\$28	\$28	\$28
P900	\$11	\$10.85	\$10.75	\$10.50	\$10.25
P910	\$11	\$10.85	\$10.75	\$10.50	\$10.25
P1000	\$35	\$34	\$33	\$33	\$33



Price List

Effective March 2013

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
R101	\$85	\$83	\$81	\$79	\$79
R110	\$85	\$83	\$81	\$79	\$79
R200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
R201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
R401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
R505	\$55	\$53	\$51	\$50	\$49
R510	\$55	\$53	\$51	\$50	\$49
R515	\$55	\$53	\$51	\$50	\$49
R520	\$55	\$53	\$51	\$50	\$49
R525	\$55	\$53	\$51	\$50	\$49
R535	\$55	\$53	\$51	\$50	\$49
R700	\$35	\$34	\$33	\$33	\$33
R750	\$35	\$34	\$33	\$33	\$33
R800	\$30	\$29	\$28	\$28	\$28
R801	\$30	\$29	\$28	\$28	\$28
R900	\$11	\$10.85	\$10.75	\$10.50	\$10.25
R910	\$11	\$10.85	\$10.75	\$10.50	\$10.25
R1000	\$35	\$34	\$33	\$33	\$33
S101	\$105	\$103	\$101	\$99	\$99
S110	\$105	\$103	\$101	\$99	\$99
S505	\$55	\$53	\$51	\$50	\$49
S510	\$55	\$53	\$51	\$50	\$49
S515	\$55	\$53	\$51	\$50	\$49
S535	\$55	\$53	\$51	\$50	\$49
S700	\$35	\$34	\$33	\$33	\$33
S750	\$35	\$34	\$33	\$33	\$33
S800	\$30	\$29	\$28	\$28	\$28
S801	\$30	\$29	\$28	\$28	\$28
S910	\$12	\$11.85	\$11.75	\$11.50	\$11.25
S1000	\$35	\$34	\$33	\$33	\$33
T5100	\$50	\$49	\$48	\$48	\$48
T5105	\$50	\$49	\$48	\$48	\$48
T5120	\$50	\$49	\$48	\$48	\$48
T5125	\$50	\$49	\$48	\$48	\$48
T5140	\$50	\$49	\$48	\$48	\$48
T5145	\$50	\$49	\$48	\$48	\$48
T5160	\$50	\$49	\$48	\$48	\$48
T5170	\$50	\$49	\$48	\$48	\$48
T5175	\$50	\$49	\$48	\$48	\$48
T5180	\$50	\$49	\$48	\$48	\$48
T5185	\$50	\$49	\$48	\$48	\$48
T5190	\$25	\$24	\$24	\$24	\$24
T5220	\$7	\$6	\$5.50	\$5	\$5
T5240	\$5	\$4.60	\$4.20	\$4	\$4
T5320-6	\$12	\$11	\$10	\$10	\$10
T5340-6	\$14	\$13	\$12	\$12	\$12

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
T5400	\$35	\$34	\$33	\$33	\$33
T5405	\$35	\$34	\$33	\$33	\$33
T5410	\$35	\$34	\$33	\$33	\$33
T5415	\$35	\$34	\$33	\$33	\$33
T5420	\$35	\$34	\$33	\$33	\$33
T5425	\$35	\$34	\$33	\$33	\$33
T5430	\$35	\$34	\$33	\$33	\$33
T5435	\$35	\$34	\$33	\$33	\$33
T5440	\$35	\$34	\$33	\$33	\$33
T5445	\$35	\$34	\$33	\$33	\$33
T5500	\$30	\$29	\$29	\$29	\$29
T5510	\$30	\$29	\$29	\$29	\$29
T5520	\$30	\$29	\$29	\$29	\$29
T5530	\$30	\$29	\$29	\$29	\$29
T5540	\$30	\$29	\$29	\$29	\$29
T5550	\$30	\$29	\$29	\$29	\$29
T5580	\$10	\$9	\$9	\$9	\$9
T5590	\$180	\$175	\$175	\$175	\$175
X100	\$100	\$98	\$96	\$94	\$92
X101	\$100	\$98	\$96	\$94	\$92
X665	\$14	\$13.50	\$13	\$12.50	\$12
X750	\$35	\$34	\$33	\$33	\$33
X800	\$35	\$34	\$33	\$33	\$33
X801	\$35	\$34	\$33	\$33	\$33
X900	\$12	\$11.75	\$11.50	\$11	\$10.50
Y100	\$105	\$103	\$101	\$99	\$97
Y101	\$105	\$103	\$101	\$99	\$97
Y300	\$15	\$14.50	\$14	\$13.50	\$13
Y301	\$15	\$14.50	\$14	\$13.50	\$13
Y400	\$80	\$78	\$76	\$74	\$72
Y401	\$80	\$78	\$76	\$74	\$72
Y500	\$60	\$58	\$56	\$55	\$54
Y550	\$60	\$58	\$56	\$55	\$54
Y555	\$60	\$58	\$56	\$55	\$54
Y560	\$60	\$58	\$56	\$55	\$54
Y665	\$14	\$13.50	\$13	\$12.50	\$12
Y670	\$14	\$13.50	\$13	\$12.50	\$12
Y700	\$35	\$34	\$33	\$33	\$33
Y750	\$35	\$34	\$33	\$33	\$33
Y800	\$35	\$34	\$33	\$33	\$33
Y801	\$35	\$34	\$33	\$33	\$33
Y900	\$15	\$14.50	\$14	\$13.50	\$13
Y950	\$15	\$14.50	\$14	\$13.50	\$13
Y955	\$15	\$14.50	\$14	\$13.50	\$13
Y960	\$15	\$14.50	\$14	\$13.50	\$13



Price List

Effective March 2013

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
Z100	\$145	\$143	\$141	\$139	\$137
Z101	\$145	\$143	\$141	\$139	\$137
Z300	\$15	\$14.50	\$14	\$13.50	\$13
Z301	\$15	\$14.50	\$14	\$13.50	\$13
Z500	\$60	\$58	\$56	\$55	\$54
Z550	\$60	\$58	\$56	\$55	\$54
Z700	\$35	\$34	\$33	\$33	\$33
Z750	\$35	\$34	\$33	\$33	\$33
Z800	\$35	\$34	\$33	\$33	\$33
Z801	\$35	\$34	\$33	\$33	\$33
Z900	\$16	\$15.50	\$15	\$14.50	\$14
Z950	\$16	\$15.50	\$15	\$14.50	\$14
Z955	\$16	\$15.50	\$15	\$14.50	\$14
AA100	\$80	\$78.50	\$77	\$77	\$77
AA101	\$80	\$78.50	\$77	\$77	\$77
AA200	\$12	\$11.50	\$9.50	\$9.25	\$9
AA201	\$12	\$11.50	\$9.50	\$9.25	\$9
AA401	\$12	\$11.50	\$9.50	\$9.25	\$9
AA505	\$50	\$48	\$46.50	\$45	\$43.50
AA510	\$50	\$48	\$46.50	\$45	\$43.50
AA515	\$50	\$48	\$46.50	\$45	\$43.50
AA520	\$50	\$48	\$46.50	\$45	\$43.50
AA525	\$50	\$48	\$46.50	\$45	\$43.50
AA535	\$50	\$48	\$46.50	\$45	\$43.50
AA660	\$14	\$13.50	\$13	\$12.50	\$12
AA680	\$15	\$14.50	\$14	\$13.50	\$13
AA700	\$35	\$34	\$33	\$33	\$33
AA750	\$35	\$34	\$33	\$33	\$33
AA800	\$30	\$29	\$28	\$28	\$28
AA801	\$30	\$29	\$28	\$28	\$28
AA910	\$14	\$13	\$11	\$10.50	\$10
AA1000	\$35	\$34	\$33	\$33	\$33
AB100	\$85	\$83.50	\$82	\$82	\$82
AB101	\$85	\$83.50	\$82	\$82	\$82
AB200	\$12	\$11.50	\$9.50	\$9.25	\$9
AB201	\$12	\$11.50	\$9.50	\$9.25	\$9
AB401	\$12	\$11.50	\$9.50	\$9.25	\$9
AB505	\$50	\$48	\$46.50	\$45	\$43.50
AB510	\$50	\$48	\$46.50	\$45	\$43.50
AB515	\$50	\$48	\$46.50	\$45	\$43.50
AB520	\$50	\$48	\$46.50	\$45	\$43.50
AB525	\$50	\$48	\$46.50	\$45	\$43.50
AB535	\$50	\$48	\$46.50	\$45	\$43.50
AB660	\$14	\$13.50	\$13	\$12.50	\$12
AB680	\$15	\$14.50	\$14	\$13.50	\$13
AB700	\$35	\$34	\$33	\$33	\$33
AB750	\$35	\$34	\$33	\$33	\$33
AB800	\$30	\$29	\$28	\$28	\$28
AB801	\$30	\$29	\$28	\$28	\$28
AB910	\$14	\$13	\$11	\$10.50	\$10
AB1000	\$35	\$34	\$33	\$33	\$33

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
AC100	\$85	\$83.50	\$82	\$82	\$82
AC101	\$85	\$83.50	\$82	\$82	\$82
AC200	\$12	\$11.50	\$9.50	\$9.25	\$9
AC201	\$12	\$11.50	\$9.50	\$9.25	\$9
AC401	\$12	\$11.50	\$9.50	\$9.25	\$9
AC505	\$50	\$48	\$46.50	\$45	\$43.50
AC510	\$50	\$48	\$46.50	\$45	\$43.50
AC515	\$50	\$48	\$46.50	\$45	\$43.50
AC520	\$50	\$48	\$46.50	\$45	\$43.50
AC525	\$50	\$48	\$46.50	\$45	\$43.50
AC700	\$35	\$34	\$33	\$33	\$33
AC750	\$35	\$34	\$33	\$33	\$33
AC800	\$30	\$29	\$28	\$28	\$28
AC801	\$30	\$29	\$28	\$28	\$28
AC910	\$14	\$13	\$11	\$10.50	\$10
AC1000	\$35	\$34	\$33	\$33	\$33
AD100	\$145	\$143.50	\$142	\$142	\$142
AD101	\$145	\$143.50	\$142	\$142	\$142
AD200	\$12	\$11.50	\$9.50	\$9.25	\$9
AD201	\$12	\$11.50	\$9.50	\$9.25	\$9
AD510	\$50	\$48	\$46.50	\$45	\$43.50
AD700	\$35	\$34	\$33	\$33	\$33
AD910	\$14	\$13	\$11	\$10.50	\$10
AD1000	\$35	\$34	\$33	\$33	\$33
AE100	\$90	\$88	\$86	\$84	\$84
AE500	\$55	\$53	\$51	\$50	\$49
AE510	\$55	\$53	\$51	\$50	\$49
AE650	\$12	\$11.75	\$11.50	\$11	\$10.50
AF100	\$100	\$98	\$96	\$94	\$94
AF500	\$55	\$53	\$51	\$50	\$49
AF510	\$55	\$53	\$51	\$50	\$49
AF900	\$12	\$11.50	\$11	\$10.50	\$10
AF910	\$12	\$11.50	\$11	\$10.50	\$10
AG500	\$55	\$53	\$51	\$50	\$49
AG510	\$55	\$53	\$51	\$50	\$49
AG910	\$12	\$11.50	\$11	\$10.50	\$10
AH100	\$110	\$108	\$106	\$104	\$104
AH510	\$55	\$53	\$51	\$50	\$49
AH900	\$12	\$11.50	\$11	\$10.50	\$10
AH910	\$12	\$11.50	\$11	\$10.50	\$10



Price List

Effective March 2013

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
AJ100	\$80	\$78.50	\$77	\$77	\$77
AJ101	\$80	\$78.50	\$77	\$77	\$77
AJ201	\$12	\$11.50	\$9.50	\$9.25	\$9
AJ505	\$50	\$48	\$46.50	\$45	\$43.50
AJ510	\$50	\$48	\$46.50	\$45	\$43.50
AJ535	\$50	\$48	\$46.50	\$45	\$43.50
AJ600	\$50	\$49	\$48	\$47	\$47
AJ601	\$50	\$49	\$48	\$47	\$47
AJ650	\$12	\$11.75	\$11.50	\$11	\$10.50
AJ651	\$12	\$11.75	\$11.50	\$11	\$10.50
AJ680	\$13	\$12.75	\$12.50	\$12	\$11.50
AJ681	\$13	\$12.75	\$12.50	\$12	\$11.50
AJ700	\$35	\$34	\$33	\$33	\$33
AJ750	\$35	\$34	\$33	\$33	\$33
AJ800	\$30	\$29	\$28	\$28	\$28
AJ801	\$30	\$29	\$28	\$28	\$28
AJ910	\$14	\$13	\$11	\$10.50	\$10
AJ1000	\$35	\$34	\$33	\$33	\$33
AK100	\$90	\$88	\$86	\$84	\$84
AK101	\$90	\$88	\$86	\$84	\$84
AK200	\$12	\$11.50	\$9.50	\$9.25	\$9
AK201	\$12	\$11.50	\$9.50	\$9.25	\$9
AK505	\$50	\$48	\$46.50	\$45	\$43.50
AK510	\$50	\$48	\$46.50	\$45	\$43.50
AK515	\$50	\$48	\$46.50	\$45	\$43.50
AK520	\$50	\$48	\$46.50	\$45	\$43.50
AK525	\$50	\$48	\$46.50	\$45	\$43.50
AK700	\$35	\$34	\$33	\$33	\$33
AK750	\$35	\$34	\$33	\$33	\$33
AK800	\$30	\$29	\$28	\$28	\$28
AK801	\$30	\$29	\$28	\$28	\$28
AK910	\$14	\$13	\$11	\$10.50	\$10
AK1000	\$35	\$34	\$33	\$33	\$33
AL100	\$100	\$98	\$96	\$94	\$94
AL101	\$100	\$98	\$96	\$94	\$94
AL200	\$12	\$11.50	\$9.50	\$9.25	\$9
AL201	\$12	\$11.50	\$9.50	\$9.25	\$9
AL505	\$50	\$48	\$46.50	\$45	\$43.50
AL510	\$50	\$48	\$46.50	\$45	\$43.50
AL515	\$50	\$48	\$46.50	\$45	\$43.50
AL520	\$50	\$48	\$46.50	\$45	\$43.50
AL525	\$50	\$48	\$46.50	\$45	\$43.50
AL700	\$35	\$34	\$33	\$33	\$33
AL750	\$35	\$34	\$33	\$33	\$33
AL800	\$30	\$29	\$28	\$28	\$28
AL801	\$30	\$29	\$28	\$28	\$28
AL910	\$14	\$13	\$11	\$10.50	\$10
AL1000	\$35	\$34	\$33	\$33	\$33

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
AM100	\$130	\$128	\$126	\$124	\$122
AM101	\$130	\$128	\$126	\$124	\$122
AM200	\$12	\$11.50	\$9.50	\$9.25	\$9
AM201	\$12	\$11.50	\$9.50	\$9.25	\$9
AM505	\$50	\$48	\$46.50	\$45	\$43.50
AM510	\$50	\$48	\$46.50	\$45	\$43.50
AM515	\$50	\$48	\$46.50	\$45	\$43.50
AM700	\$35	\$34	\$33	\$33	\$33
AM910	\$14	\$13	\$11	\$10.50	\$10
AM1000	\$35	\$34	\$33	\$33	\$33
AP101	\$80	\$78	\$76	\$74	\$74
AP110	\$80	\$78	\$76	\$74	\$74
AP200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AP201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AP401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AP505	\$55	\$53	\$51	\$50	\$49
AP510	\$55	\$53	\$51	\$50	\$49
AP515	\$55	\$53	\$51	\$50	\$49
AP520	\$55	\$53	\$51	\$50	\$49
AP525	\$55	\$53	\$51	\$50	\$49
AP535	\$55	\$53	\$51	\$50	\$49
AP700	\$35	\$34	\$33	\$33	\$33
AP750	\$35	\$34	\$33	\$33	\$33
AP800	\$30	\$29	\$28	\$28	\$28
AP801	\$30	\$29	\$28	\$28	\$28
AP900	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AP910	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AP1000	\$35	\$34	\$33	\$33	\$33
AR101	\$85	\$83	\$81	\$79	\$79
AR110	\$85	\$83	\$81	\$79	\$79
AR200	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AR201	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AR401	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AR505	\$55	\$53	\$51	\$50	\$49
AR510	\$55	\$53	\$51	\$50	\$49
AR515	\$55	\$53	\$51	\$50	\$49
AR520	\$55	\$53	\$51	\$50	\$49
AR525	\$55	\$53	\$51	\$50	\$49
AR535	\$55	\$53	\$51	\$50	\$49
AR700	\$35	\$34	\$33	\$33	\$33
AR750	\$35	\$34	\$33	\$33	\$33
AR800	\$30	\$29	\$28	\$28	\$28
AR801	\$30	\$29	\$28	\$28	\$28
AR900	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AR910	\$11	\$10.85	\$10.75	\$10.50	\$10.25
AR1000	\$35	\$34	\$33	\$33	\$33



Price List

Effective March 2013

CAT. No.	PRICE PER EACH AT QUANTITY				
	1	6	12	24	36
AS101	\$105	\$103	\$101	\$99	\$99
AS110	\$105	\$103	\$101	\$99	\$99
AS505	\$55	\$53	\$51	\$50	\$49
AS510	\$55	\$53	\$51	\$50	\$49
AS515	\$55	\$53	\$51	\$50	\$49
AS535	\$55	\$53	\$51	\$50	\$49
AS700	\$35	\$34	\$33	\$33	\$33
AS750	\$35	\$34	\$33	\$33	\$33
AS800	\$30	\$29	\$28	\$28	\$28
AS801	\$30	\$29	\$28	\$28	\$28
AS910	\$12	\$11.85	\$11.75	\$11.50	\$11.25
AS1000	\$35	\$34	\$33	\$33	\$33

SHIPPING:

1. Free 2nd day shipping for orders \$150+,
2. Free standard overnight shipping for orders \$250+
3. Shipping surcharges apply to residential & out of service areas.
4. Additional shipping charges apply to Hawaii, Alaska & International destinations.

NOTE:

Prices are subject to change without notice and are in U.S. Dollars

RETURN POLICY:

1. Call for Return Authorization Number before shipping.
2. Packages must be un-opened, un-damaged & clean.
3. Must have been purchased within the last 30 days.
4. Include a copy of invoice.
5. Do not staple or tape the packages.
6. All returns MUST be shipped in a box, not envelope or padded envelope.

[illegible]

[illegible]

Lifetime Warranty

DiamoDent offers “Lifetime Warranty” on all restorative components. DiamoDent warrants that its products shall be free from defects in material and/or workmanship. This warranty applies only to the original purchaser. DiamoDent will at its option, either repair, replace or issue credit for such products. Purchaser assumes all risks and liability resulting from the use of these products. DiamoDent makes no other warranties, express or implied, except as set forth above.

In the event of a product defect, please immediately notify DiamoDent’s Customer Service Department prior to returning the product. Product must be sterilized prior to return.

DiamoDent®

1577 North Harmony Circle
Anaheim, CA 92807 U.S.A.

WWW.DIAMODENT.COM

Customer Service
Monday – Friday 7am-5pm, Pacific Time
888-281-8850
714-701-0171
714-701-0176 Fax

All products are approved by FDA

All products are manufactured in U.S and Europe

Prices are subject to change without notice

Trademarks

DiamoDent is registered trademark of DiamoDent Inc.

ComScrew & TiteFit are trademarks of DiamoDent Inc.

BIOMET, Osseotite XP, Miniplant, MicroMiniplant, Certain & 3i are registered trademarks of Implant Innovations Inc.

NobelReplace, Tapered & Straight Groovy, NobelSpeedy Replace, Replace Select Tapered & Straight, Branemark,

Steri-Oss, UniGrip, WP, RP & HL are registered trademarks of Nobel Biocare AB

Lifecore, Keystone, RENOVA, RESTORE, SD, RD, WD, SDI & RDI are registered trademarks of Keystone Dental.

BioHorizon is registered trademark of BioHorizons USA

Screw-Vent & Zimmer are registered trademarks of Zimmer

Frialit-2, XiVE & Friadent are registered trademarks of Friadent

Solid Abutment, SCS, ITI & Straumann are registered trademarks of Straumann

RePlant, RePluCs, ReActive, Legacy & Implant Direct are

registered trademarks of Implant Direct LLC

L1004 Rev C

